

Adult Continuous Suction Regulators

In Suction Management, We Are The Calm Inside the Storm.

Trusted Performance Across All Your Departments.

- Intensive Care Units
- Emergency Department

Features

- Factory 12-year warranty
- Low demonstrated cost of ownership
- Patented Self-Clearing™ technology ensures years of trouble-free operation
- No calibration required



Linear Gauge

- As suction increases, the piston goes down indicating the level of suction with a 180° field of view
- Color coding readily shows common levels of suction

Adjusting Knob

- Turn COUNTER-CLOCKWISE to increase the suction setting
- Turn CLOCKWISE to decrease suction setting
- Similar function to flowmeters

Mode Selector Valve

- Easily switch between operating modes OFF-REG-LINE
- Regulated Suction range is 20 - 200 mmHg



Description	Part Number	Ordering
3800 Suction Regulator	3800	1 unit