

OR Suction Regulators

In Suction Management, We Are
The Calm Inside the Storm.

Dependable Clinical Use:

- Operating Room
- Procedure Rooms
- High-Vacuum, High Flow Suction
- Advanced interior-regulating mechanism prevents flooding and contamination - the two most common causes for regulator failure
- Eliminates the need for bulky, costly portable suction units
- sterile reprocessing
- Allows for easy-install trap bottle and still be attached to roll stand

Distinctive Features

- High-Flow / >150lpm with fittings
- Factory 12-year warranty
- Patented Self-Clearing™ technology ensures years of trouble-free operation
- Sterilizable - VSS270™



Linear Gauge

- As suction increases, the piston goes down indicating the level of suction with a 180 field of view
- Color coding readily shows common levels of suction

Adjusting Knob

- Turn COUNTER-CLOCKWISE to increase the suction setting
- Turn CLOCKWISE to decrease suction setting
- Similar function to flowmeters

Mode Selector Valve

- Easily switch between operating modes OFF-REG-LINE
- Regulated Suction range is 76 - 760 mmHg



Description	Part Number	Ordering
3840 OR Suction Regulator	3840	1 unit