

Care4TM Suction Regulators

Suction Regulators from

BOEHRINGER[®]

**In Suction Management, We Are
The Calm Inside the Storm.**

**The Only 4-mode Suction Regulator
for All Your Clinical Suction Needs.**

- Critical Care Units
- Emergency departments
- Operating rooms
- Recovery rooms
- Step-down units
- Procedure areas

Features

- Fully-regulated continuous and intermittent suction
- High vacuum, safely accessible via integrated locking mechanism requiring clinician activation
- Full 12-year factory warranty
- Self-ClearingTM technology
- Sterilizable – VSS270



Linear Gauge

- As suction increases, the piston goes down indicating, the level of suction with a 180° field of view
- Color coding readily shows common levels of suction

Adjusting Knob

- Turn COUNTER-CLOCKWISE to increase the suction setting
- Turn CLOCKWISE to decrease suction setting
- Similar function to flowmeters



Mode Selector Valve

- Easily switch between operating modes OFF-INTERMIT-CONT-LINE
- LINE suction mode is engaged by pulling the control knob out and rotating to LINE position
- Pull outwards and towards you for LINE
- This prevents accidental activation of the LINE function



▶ SCAN FOR VIDEO

Description	Part Number	Ordering
Care4 TM Suction Regulator	3844	1 unit