

# VAVD Controller & VacPac®

Suction Regulators from


**BOEHRINGER**®

In Suction Management, We Are  
The Calm Inside the Storm.

Relied On By 600+ Heart Programs Across the US

## Features

- Excess Negative Pressure Relief
- Positive Pressure Relief
- Regulated Range Conforms to Modern VAVD Standards
- Integrated Waste Anesthetic Gas Disposal
- Dedicated Vacuum Assist Protection Circuit (VacPac®)



PART #	DESCRIPTION	ORDERING
3931	VAV Controller with WAGD and VacPac®	1 Unit
3948	VacPac® Disposable	1 Box (10 Sets)
3949	VacPac® Disposable	1 Box (40 Sets)