

# MRI Suction Regulator

In Suction Management, We Are  
The Calm Inside the Storm.

## Clinical Use Area:

- MRI Suite

## Distinctive Features

- Flood- and shock-resistant gauge needs no calibration
- High flow capacity
- Industrial strength diaphragm
- Factory qualification testing available
- 12-Year warranty
- Sterilizable – VSS270™



### Linear Gauge

- As suction increases, the piston goes down indicating, the level of suction with a 180 field of view
- Color coding readily shows common levels of suction

### Adjusting Knob

- Turn COUNTER-CLOCKWISE to increase the suction setting
- Turn CLOCKWISE to decrease suction setting
- Similar function to flowmeters

### Mode Selector Valve

- Integrated locking mechanism requires direct clinician intervention to engage full-line suction- preventing accidental activation of "line" function
- Regulated Suction range is 20-200 mmHg
- Modes: OFF-REG-LINE



Description	Part Number	Ordering
MRI Suction Regulator	3708	1 unit