

# Dual Suction Regulator Mount

In Suction Management, We Are  
The Calm Inside the Storm.



**EXPAND YOUR ACCESS TO SUCTION:**  
Our Dual Suction Regulator Mount allows for two suction regulators to be used simultaneously from a single vacuum outlet for interventions such as:

- In-line suction for airway maintenance
- Nasogastric suction to lessen the chance of aspiration
- Oral Suction to establish an airway for oral care
- Continuous suction for external female catheter use
- Full-line suction for emergency airway clearance, bedside procedures, etc.

Description	Part Number	Ordering
Dual Suction Regulator Mount	9063T	1 unit